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
INTERCEPT PLATELET PROCESSING SET with DUAL STORAGE CONTAINERS

For use with INTERCEPT Illuminator

Each set is wrapped in a tamper-evident package and includes one 17.5mL 3mM amotosalen hydrochloride solution container (Formula : Amotosalen HCl 101mg - Natr. chlorid. 924mg - Aqua ad iniect. ad 100mL), one illumination container, one container with Compound Adsorption Device (CAD), two INTERCEPT platelet storage containers. The set is sterilized by a combination of steam and radiation.

Indications and Usage

Intended Use

The set is used with an INTERCEPT Illuminator to inactivate a broad spectrum of viruses, bacteria, and parasites as well as contaminating donor leukocytes in platelet components. This process for treatment of platelet components is intended to reduce the risk of transfusion-associated transmission of viruses, bacteria, and parasites, and the risk of adverse effects due to transfusion of contaminating donor leukocytes.

Indications for Use

INTERCEPT Platelets are indicated for transfusion support of patients requiring platelet transfusions according to clinical practice guidelines. Any type of thrombocytopenia or qualitative disorder resulting from disease, therapy, or injury can be supported with INTERCEPT Platelets. INTERCEPT treatment may be used as an alternative to gamma irradiation for prevention of transfusion-associated graft-versus-host disease (TA-GVHD). INTERCEPT treatment may be used in place of CMV testing and leukoreduction for prevention of transfusion transmitted CMV infection. INTERCEPT Platelets are not clinically different from untreated platelets and are infused according to standard platelet infusion methods.

INTERCEPT platelets suspended in additive solution or in 100% plasma may be stored up to 7 days from time of collection. Treated platelets must be stored at 20-24°C with continuous agitation. Any extension of platelet storage time should be evaluated according to applicable local policies and regulations.

Platelet additive solutions approved for use with INTERCEPT: InterSol, SSP+, T-PAS+, Grifols PAS III M.

Contraindications

Use of INTERCEPT Platelets is contraindicated in patients with a history of allergic response to amotosalen or psoralens.

Precautions

Do not use if : tamper-evident package has been opened; signs of deterioration are visible; fluid path closures are loose or not intact; cannulae are broken or there is no fluid in amotosalen solution container. Do not store above 25°C. Do not vent. Do not freeze. Protect the pack and tubing from sharp objects. Unused sets in open aluminium foil may be kept 20 days at room temperature by folding and securing open end of aluminium foil. Units removed from the aluminium foil must be used within 8 hours. Keep set in light-protective package until time of use. Protect from direct sunlight and strong UVA light source.

Set is single use only. Do not reuse. Do not resterilize. This product is not designed for reuse. Misuse can result in adverse reactions, including severe illness and possibly death.

All the following conditions must be met for pathogen inactivation:

- Platelets must be prepared with additive solution for platelet re-suspension in the volume range specified in **Table 1**.
- Platelet content, volume, and red blood cell (RBC) count must be within the ranges specified in **Table 1**.
- Platelets mixed with amotosalen must be exposed to UVA light dose from the INTERCEPT Illuminator. No other source of UVA light may be used.
- Platelets collected on Day 0 must be exposed to UVA light by end of Day 1.
- After illumination, platelets must be agitated in CAD container in accordance with the durations specified in **Table 2**.

This process is designed to be a closed system. Treatment with INTERCEPT Blood System does not replace applicable standards for processing in open and closed systems. If there is a leak in the set during processing, platelet product must be discarded.

Warnings : Amotosalen in contact with skin may result in photosensitisation in the presence of ultraviolet light. If skin exposure occurs, flush exposed skin copiously with water. Sterile connecting device (SCD) should be used according to manufacturer's instructions for use.

Notes to Physicians

While laboratory studies of amotosalen processing with UVA light have shown a reduction in levels of certain viruses, bacteria and parasites, there is no pathogen inactivation process that has been shown to eliminate all pathogens.

INTERCEPT platelets components should not be prescribed to neonatal patients treated with phototherapy devices that emit a peak energy wavelength less than 425 nm, and/or have a lower bound of the emission bandwidth <375 nm, due to the risk of erythema resulting from potential interaction between ultraviolet light (below 400 nm) and residual amotosalen.

Instructions for Use

Materials Needed: One (1) INTERCEPT Platelet Processing Set with Dual Storage Containers.

Equipment Needed : INTERCEPT Illuminator, Sterile Connecting Device (SCD), Tube Sealer, Flatbed Agitator.

Temperatures should be controlled to meet applicable regulations for platelet processing.

Used and unused INTERCEPT Processing sets should be discarded like any used blood containers, as biohazardous waste.

Process Steps

A- Preparation of Platelets

Platelets can be prepared in additive solution or in 100% plasma (**Table 1**). Platelet products within the following ranges have been shown to be acceptable for use with this processing set.

Table 1

Suspension Medium		Platelet Count	Volume	RBC Content
Plasma Content	Additive Solution Content			
32 - 47%	53 - 68%	2.5 - 7.0 x 10 ¹¹	300 - 420 mL	<4 x 10 ⁶ RBC/mL
		7.1 - 8.0 x 10 ¹¹	375 - 420 mL	
100%	0%	2.5 - 7.0 x 10 ¹¹	300 - 420 mL	<4 x 10 ⁶ RBC/mL
		7.1 - 8.0 x 10 ¹¹	375 - 420 mL	

B- Amotosalen Addition to Platelets

1. Remove set from package. Unwrap only illumination container from organizer.
2. Weld tubing from platelet container to amotosalen container tubing using SCD.
3. Label set containers using appropriate donation identification. See warning in Section C. While labeling storage containers, separate them to ensure they do not adhere to one another.
4. Hang platelets and break both cannulae on amotosalen container.
5. Allow platelets and amotosalen solution to flow into illumination container, marked by number "1".
6. Ensure that platelets are completely transferred to illumination container by expressing air from platelet and amotosalen container into illumination container.
7. When platelets are in illumination container, mix thoroughly by gentle agitation to ensure complete mixing of amotosalen and platelets.
8. Express air from platelets into amotosalen container.
9. Express a small amount of platelet and amotosalen mixture into tubing, filling about 4cm of tubing.
10. Seal tubing between illumination container and amotosalen container, so that tubing is no longer than approximately 4cm from illumination container inlet port.
Warning : During illumination, tubing must be held within large compartment of illumination tray.
11. Remove and discard empty platelet and amotosalen containers.

C- Illumination

Illuminate platelets.

Refer to INTERCEPT Illuminator Operator's Manual for complete illumination instructions for use.

Warning : Platelets in illumination container must be within large compartment of illuminator tray in order for inactivation to occur. The process requires unimpeded light transmission through tray and illumination container with platelets. No labels or other material should be on this area. Tray must be clean. Labels should be placed on illumination container flap only. Illumination container should lay flat.

D- Agitation with CAD

Warning : Do not fold or bend CAD.

1. Unwrap set from organizer.
2. Hang platelets, break cannula and allow platelets to flow into CAD container, marked by number "2".
3. Express air from CAD container into illumination container.
4. Seal tubing close to inlet port of CAD container.
5. Remove and discard empty illumination container.
6. Place CAD container on flatbed agitator for the duration specified in **Table 2**.

Table 2

Suspension Medium	CAD Agitation Duration
Additive Solution (53 - 68%)	6 to 16 hours
Plasma (100%)	16 to 24 hours

E- Transfer to Storage Containers

1. Remove platelets from agitator and hang platelets.
2. Close clamp on platelet sampling pouch.
3. Break cannula and allow platelets to flow into each of the two storage containers, marked by number "3."
4. Express air from storage containers into CAD container and clamp above the Y fitting.
5. Redistribute volume between the storage containers as needed. Note: Steps 4 and 5 may be performed in the reverse order.
6. Ensure appropriate donor identification is attached to each filled platelet storage container.
7. Seal tubing to each storage container close to outlet of the Y fitting.
For platelets in 100% plasma: If the input platelet count to the INTERCEPT process is greater than or equal to 5.3 x 10¹¹, the product should be split between two containers for storage.
8. Remove and discard empty CAD container. Place platelets on platelet agitator.
The INTERCEPT Platelet process is now complete.

F- Sampling Platelet Product (optional)

1. Mix INTERCEPT platelets by gently agitating platelet storage container.
2. Open clamp to platelet sampling pouch and squeeze several times.
3. Allow sampling pouch to fill with platelets. Seal tubing.
4. Remove sample pouch.
5. Transfer sample to appropriate laboratory tube immediately.

Caution: DEHP is known to be released from polyvinyl chloride (PVC) medical devices; increased leaching can occur with extended storage or increased surface area contact. The INTERCEPT processing sets only have tubing components, container ports and, if included, an in-line filter that contain PVC; all containers and other parts are PVC-free. During use of this processing set, blood components are in contact with PVC for a brief period of time (approx. <15 minutes). Based on limited surface area contact and minimal contact time, DEHP levels in blood components after use of the processing set are estimated to be well below those resulting from other medical applications containing PVC tubing (e.g. hemodialysis, intravenous fluid administration, extracorporeal membrane oxygenation and cardiopulmonary bypass procedures). The risks associated with DEHP released to the blood components must be weighed against the benefits of therapeutic transfusion and inactivation of harmful viruses, bacteria and other pathogens.

Cerus, INTERCEPT, and INTERCEPT Blood System are trademarks of Cerus Corporation

InterSol is a trademark of Fenwal, Inc.

SSP+ is a trademark of MacoPharma

T-PAS+ is a trademark of Terumo BCT, Inc.

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Grifols PAS III M is a trademark of Grifols