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
INTERCEPT PROCESSING SET for PLASMA

For use with INTERCEPT Illuminator

Each set is wrapped in a tamper-evident package and includes one 15mL 6mM amotosalen hydrochloride solution container (Formula : Amotosalen HCl 203mg - Natr. chlorid. 924mg - Aqua ad inject. ad 100mL), one illumination container, one compound adsorption device (CAD), three INTERCEPT Plasma storage containers. The set is sterilized by a combination of steam and radiation.

Indications and Usage

Intended Use:

This set is used with an INTERCEPT Illuminator to inactivate a broad spectrum of viruses, bacteria, parasites as well as contaminating donor leukocytes in plasma.

Indications for Use:

INTERCEPT plasma is indicated for support of patients requiring plasma transfusions or therapeutic plasma exchange, according to clinical practice guidelines. Clinical trials in patients have demonstrated that plasma treated with the INTERCEPT Blood System was well tolerated and retained therapeutic efficacy comparable to conventional fresh frozen plasma. INTERCEPT Plasma may be used to treat single coagulation factor or antithrombotic protein deficiencies for which no concentrates are available, as well as multiple coagulation factor and antithrombotic protein deficiencies. INTERCEPT Plasma may also be used for plasma exchange for thrombotic thrombocytopenic purpura (TTP). INTERCEPT treatment may be used as an alternative to gamma irradiation for prevention of transfusion-associated graft-versus host disease (TA-GVHD). INTERCEPT treatment may be used in place of CMV testing and leukoreduction for prevention of transfusion-transmitted CMV infection.

INTERCEPT plasma may be stored from the time of collection for 12 months between -18°C and -25°C or for 24 months below -25°C, in compliance with applicable procedures and regulations.

Plasma photochemically treated with the INTERCEPT Blood System may be stored and transfused according to standard methods for frozen plasma. Thawed INTERCEPT Plasma that has been stored at 2-6°C can be used up to 5 days. As with all plasma products, clinical use should consider that labile coagulation factors decline during post thaw storage.

Contraindications

Use of INTERCEPT Plasma is contraindicated in patients with a history of allergic response to amotosalen or psoralens.

Precautions

Do not use if: tamper-evident package has been opened; signs of deterioration are visible; fluid path closures are loose or not intact; cannulae are broken or there is no fluid in amotosalen solution container.

Do not store above 25°C. Do not vent. Do not freeze. Protect the pack and tubing from sharp objects.

Unused sets in open aluminium foil may be kept 20 days at room temperature by folding and securing open end of aluminium foil. Units removed from the aluminium foil must be used within 24 hours.

Keep set in light-protective package until time of use. Protect from direct sunlight and strong UVA light source.

Set is single use only. Do not reuse. Do not sterilize. This product is not designed for reuse. Misuse can result in adverse reactions, including severe illness and possibly death.

All the following conditions must be met for pathogen inactivation:

- Plasma volume and red blood cell (RBC) content must be within the range specified in Table 1.
- Plasma mixed with amotosalen must be exposed to UVA light dose from INTERCEPT Illuminator. No other source of UVA light may be used.
- Plasma must be passed through the CAD by gravity flow process after illumination.
- For fresh plasma, the entire process from collection to freezing must be completed within a timeframe of 20 hours or according to local regulations.
- Fresh frozen plasma that has been thawed in accordance with local regulations must be immediately treated with the INTERCEPT Blood System for plasma and transfused or refrozen promptly after treatment.

This process is designed to be a closed system. Treatment with INTERCEPT Blood System does not replace applicable standards for processing in open and closed systems. If there is a leak in the set during processing, plasma product must be discarded.

Warnings : Amotosalen in contact with skin may result in photosensitisation in the presence of ultraviolet light. If skin exposure occurs, flush exposed skin copiously with water. Sterile connecting device (SCD) should be used according to manufacturer's instructions for use.

Notes to Physicians

While laboratory studies of amotosalen processing with UVA light have shown a reduction in levels of certain viruses, bacteria and parasites; there is no pathogen inactivation process that has been shown to eliminate all pathogens.

INTERCEPT plasma components should not be prescribed to neonatal patients treated with phototherapy devices that emit a peak energy wavelength less than 425 nm, and/or have a lower bound of the emission bandwidth <375 nm, due to the risk of erythema resulting from potential interaction between ultraviolet light (below 400 nm) and residual amotosalen.

Instructions for Use

Materials Needed: One (1) INTERCEPT Processing Set for Plasma.

Equipment Needed: INTERCEPT Illuminator, Sterile Connecting Device (SCD), Tube Sealer.

Temperatures should be controlled to meet applicable regulations for plasma processing. Used and unused INTERCEPT sets should be discarded like any used blood containers, as biohazardous waste.

Process Steps

A-Plasma for Processing with the INTERCEPT Blood System

Plasma products within the following ranges have been shown to be acceptable for use with this processing set.

Table 1

Volume	RBC content
385-650 mL	<4 × 10 ⁶ RBC/mL

B- Amotosalen Addition

1. Remove set from package.
2. Weld tubing from plasma container to amotosalen container tubing using SCD.
3. Disassemble set from organizer and remove rubber band.
4. If two plasma units will be produced by the INTERCEPT process, heat seal and remove one storage container.
5. Label set containers using appropriate donation identification. See warning in section C. While labeling storage containers, separate them to ensure they do not adhere to one another.
6. Hang plasma container, ensuring that set containers/components do not come in contact with floor. Break both cannulae on amotosalen container.
7. Allow plasma and amotosalen solution to flow into illumination container marked by number "1".
8. Ensure that plasma drains completely from initial plasma container into illumination container.
9. Express air from the illumination container into the amotosalen container.
10. When air is removed and plasma has fully drained back into the illumination container, mix illumination container thoroughly by gentle agitation to ensure complete mixing of amotosalen and plasma.
11. Express a small amount of plasma and amotosalen mixture into tubing, filling at least 4cm of tubing.
12. Seal tubing between illumination container and amotosalen container within the 4cm. **Warning :** During illumination, tubing must be held within large compartment of illumination tray.
13. Remove and discard initial plasma container, amotosalen container and excess tubing.

C- Illumination

Illuminate plasma.

Refer to INTERCEPT Illuminator Operator's Manual for complete illumination instructions for use. **Warning :** All plasma, both in illumination container and tubing, must be within large compartment of illuminator tray in order for inactivation to occur. The process requires unimpeded light transmission through tray and illumination container with plasma. No labels or other material should be on this area. Tray must be clean. Labels should be placed on illumination container flap only. Illumination container should lay flat.

D- Processing with Compound Adsorption Device (CAD)

1. Hang illumination container, allowing CAD to hang freely, with storage containers kept in an inverted position.
2. Close clamp on bypass line; ensure all other clamps are open.
3. Break cannula on illumination container and allow plasma to flow through CAD marked with the number "2" into storage containers marked with the number "3".
4. Once plasma has emptied from illumination container and passed through CAD, close clamp on tubing leading from the CAD and open clamp on the bypass line.
5. Hang CAD together with illumination container.
6. Remove storage containers from tab on CAD and allow them to hang ports up.
7. Completely express air from storage containers through bypass line.
8. Close clamp on the bypass line and open the clamps on storage containers.
9. Open clamp below the CAD, allowing plasma to drain into storage containers.
10. Ensure that the storage container tubing contains plasma and no air. Close clamps on storage containers.
11. Re-distribute plasma volume between storage containers, if necessary.
12. Ensure appropriate donor identification is attached to each filled plasma storage container.
13. Disconnect storage containers from set by heat-sealing, allowing sufficient tubing length for segments.
14. Discard CAD and illumination container. The INTERCEPT Plasma process is now complete.
15. Seal tubing as appropriate for making segments as needed.
16. Follow internal procedures for freezing plasma.

Caution: DEHP, when used as a plasticizer, is known to be released from polyvinyl chloride (PVC) medical devices; increased leaching can occur with extended storage or increased surface area contact. The INTERCEPT processing sets only have the cannulae and Y-junction connectors that contain PVC with DEHP; all containers and other parts are DEHP-free. During use of this processing set, blood components are in contact with DEHP for a brief period of time (approx. <15 minutes). Based on limited surface area contact and minimal contact time, DEHP levels in blood components after use of the processing set are estimated to be well below those resulting from other medical applications containing PVC with DEHP tubing (e.g. hemodialysis, intravenous fluid administration, extracorporeal membrane oxygenation and cardiopulmonary bypass procedures). The risks associated with DEHP released to the blood components must be weighed against the benefits of therapeutic transfusion and inactivation of harmful viruses, bacteria and other pathogens.

Cerus, INTERCEPT, and INTERCEPT Blood System are trademarks of Cerus Corporation

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„Cerus“, INTERCEPT ir „INTERCEPT Blood System“ yra „Cerus Corporation“ prekių ženklai.